

# P/S SERIES PRINTING THROUGH PC GUIDE

- Patient monitor needs to be connected to the computer through network (Example of IP address: Monitor – 192.168.1.XXX; PC – 192.168.1.XXX). To test whether the connection is successful, you can check it by ping it in the terminal.
- 2. Ensure PC is connected to the printer (Printer driver needs to be installed).

### Instructions:

- 1. Installation package will be provided, extracts the .Zip file into any path.
- 2. After extraction, open the <u>print server v1.1</u> folder, open <u>exe</u> folder and look for <u>miniprogram.exe</u>. Double–click on <u>miniprogram.exe</u> to run the program.
- 3. The program only supports single instances, and if try multiple run, a warning will pop up to alert you that's already running, operating process can be viewed in the taskbar. Warnings as shown in the following figure:

🔳 Warn	ing ×
	Biolight Background Print Program has run !
	ОК

If you don't see it, it's probably hidden. Click the Show Hidden Tasks button in the taskbar. Please ensure this program is running when start the PC.

### Printer settings:

4. Printer needs to be configured before working properly. Upon first operation of the program, printer setup window will pop up and require user to select a printer from list. As shown in the following figure:

🖨 Print Settings	—		$\times$		
Please select a valid p	rinter from	the li	st :		
Canon LBP2900Plus			•		Click to select printer
Please select a print s	tyle :				
One-side Print			•	<b></b>	Click to set printing style
					* Default: Single-sided

5. If there are multiple printers installed, you can try to print any pdf (such as test.pdf) to clarify the working printer, the software does not specify the working printer.



# Monitor settings:

- 1. Printing program act as a server, and IP address of PC need to be set on patient monitor.
- After obtaining the IP of the PC, go to Maintenance of patient monitor and input Maintenance password (785623) → Printing → Printer, set the printing service IP, as shown in the following figure:

	Maintenance										
Device Location	Patient Management	Alarm	CAA	Module	Print	Unit	Othe	•			
Recorder	Print Server IP						19	2.168	.0.23		
Printer											

3. After setting up, **Freeze waveform** or **Real-time waveform** can be printed from the P Series Monitor by pressing respective tab as shown in the following figure:







# Follow the step to set up report.

👔 4 Jimmy	y Zhao									Default 💢	💥 🦛	🕂 🚅	💭 04:05:42 PM	
Pediatric ECG	Leads Off								🖄 Alarm	Pause 3s				
II Diagnose	12.5mm/s Notch 5							ECC	3 Leads Off					
xt														
1mV											_		SOI : 0	
				Main Menu	ı		E E E E E E E E E E E E E E E E E E E							
	Disp	lay	Param	eter	Configura	ation	Report						"II	
, 	Screen	Select	Parameter	Switch	Load Curren	t Config	Report Set	tup						
1mV	Screen	Layout	Setu	p	Set Default	Config	Record Set	tup		38.0 36.0	38. 36.	0	ΔΤ	
RESP x1 RA	R <sup>®</sup> Quick Keys Param Color		Color	Import User Config					RESP					
	Oth	er			Export User	Config					_	•	ECG	
Dist	Privacy Mode				Save User (	Config								
Pleth	Night	Mode			Delete User Config					Sp02 ×			PI	
					New Patient	t Config							bpm	
NIBP					Boot Co	nfig			E	NEWS		=		
mmHg 120				1		1			v		? )	Inp	ut Data	
70				0.	0				м	@				
		22	2		<u>م</u> ۲	₽	46	9, <del>9</del>	<u>B</u>	Q.2				
	Review	Alarm Reset	Alarm Pause	Screen Setup	Alarm Setup	Voice Assistant	Intubation Status	Venipuncture	Patient Managem	ent NIBP Start/S	top		Main Menu	

4 Jimmy Pediatric Batt	<sup>y Zhao</sup> t <b>ery Not Exis</b> t	t								Default 🔀	ţ	u 💭	04:05:53 PM
II Diagnose	12.5mm/s Notch 5	DHz						ECG Les	ads Off EC (	168			
← ●				Report	Setun				×				PVCs :-
Realtime Report				¢	Auto				T1				" _==100
				e	5.25 mm/s	Sele	ect Waveform		э <sup>г</sup>				
Speed				25 mm/s	12.5 mm/s	orms				38.0 36.0			
				20 1111/0	25 mm/s	form							
				5	50 mm/s	Sel	ect Waveform		RES	;P			
						36				_		✿ ECG	
									2				
									SpC				
									*			E PR	
									2	2		2	
									NE	ws 2			
	_	_	_	_	_	_	_	_		:		Input Data	
									Print	@:			
		Z)	\$		<b>A</b> 1	Ð	<b>ب</b>	₹. <sup>6</sup>	Patient	et.			
	Review	Alarm Reset	Alarm Pause	Screen Setup	Alarm Setup	Voice Assistant	Intubation Status	Venipuncture	Management	NIBP Start/Stop			Main Menu

Waveform speed can be set as: Automatic, 6.25, 12.5, 25, 50mm/s.



Pediatric ECG	y Zhao Leads Off 12.5mm/s X Iteleh 5	DH:		Report	Setup			ECGLe	ads Off E	Default Default	• 🕂 🖟	2	04:05:48 PM sqi: 0 pvcs :-
Realtime Report									T				" <b>1</b>
						Sel	ect Waveform						
Speed				25 mm/s ▹	O Current Wavefo	orms				38.0 <b></b>	38.0 36.0		
					Selected Wave	form				FSP			
				(		Sel	lect Waveform			m			
												ECG	
									S 76	p02			
										~		PR bpm	
									Ν	IEWS			
												Input Data	
									Print	@:-			
		×1	~		<b>∆</b> ſ	Ð	ት	8.º	Batiant	Q.,			
	Review	Alarm Reset	Alarm Pause	Screen Setup	Alarm Setup	Voice Assistant	Intubation Status	Venipuncture	Management	NIBP Start/Stop			Main Menu

Printing waveforms can be selected base on users preference.

4 Jimmy Zhao Pediatric ECG Leads Off	h 50Hz						ECG Leads Off	ECG bpm	Default 💥	æ	04:06:57 PM
<ul><li>← ●</li></ul>			Report S	Setup				×			PVCs :
Realtime Report								т1			
		Choices			Selected	reform		J.			
	Pleth								8.0	38.0	ΔΤ
Speed	PlethL		Add	aVR							
	RESP		< Dele	aVF				RES	P		
	C02		< Delete	All V2		veform		rpm			
	ICG			V3					_		O ECG
	BIS		Move To	Top Vb							
	AA		Move To	Up V5				SpO	2		
	02		Move To I	Down V6				×			PI
	N20		Move To B	ottom							PR bpm
	PAW		WOVE TO D					10			
	FLOW		Default Se	aVL				NE	NIS		
	VOL			Va					ີ່ <u>ງ</u>		
											Input Data
								Print	@;		
	20	~		<u>م</u> ړ	Ą	46	9.5°	Ê.	64.5		
Review	Alarm Reset	Alarm Pause	Screen Setup	Alarm Setup	Voice Assistant	Intubation Status	Venipuncture	Patient Nanagement	NIBP Start/Stop		Main Menu



## Freeze 12 leads and print



#### 2/2

# Standard screen freeze and print





#### **Real time print**





Trend table print



Patient ID : A1234567 Age : 24years Name : JimmyBed : ICU BAY 2Height / Weight : 120 cm / 45 kg Gender : Male

Birth Date : 1998-06-19 Paced : No

2022-09-23	19:28:20	19:29:20	19:30:20	19:31:20	19:32:20	19:33:20	19:34:20	19:35:20	19:36:20
bpm	60	60	60	60	60	60	60	60	60
Sp02 %	99	99	99	99	99	99	99	99	99
PI %	7.45	7.45	7.45	7.45	7.45	7.45	7.45	7.45	7.45
SpO2L-SpO2 %									
SpO2L-PI %									
SpO2L-SpCO %									
SpO2L-SpOC ml/dl									
SpO2L-SpMet %									
Sp02L-PVI %									
SpO2L-SpHb g/dL									
RR bpm	20 Source: CO2	20 Source: CD2	20 Source: CO2						
NIBP mmHg	/() @-∵-:-	/() @∹∹∹	/() @∹∹∹	/() @∹∹∹	/() @∹∹∹	/() @-∵	125/84(96) @19:33:25	/() @∹∹∹	/() @∹∹∹
PR bpm	60 Source: SpO2	60 Source: Sp02	60 Source: SpO2						
P1 mmHg	122/81(99)	122/81(99)	122/81(99)	122/81(99)	122/81(99)	122/81(99)	122/81(99)	122/81(99)	122/81(99)
P2 mmHg	()	()	()	()	()	()	()	()	()
ART mmHg	123/82(99)	123/82(99)	123/82(99)	123/82(99)	123/82(99)	123/82(99)	123/82(99)	123/82(99)	123/82(99)
CVP cmH20	()	()	()	()	()	()	()	()	()
PA mmHg	()	()	()	()	()	()	()	()	()
PAWP mmHg	 @-	 @	 @-	 @-	 @-	@-	 @-	 @-	 @-
RAP mmHg	()	()	()	()	()	()	()	()	()
LAP mmHg	()	()	()	()	()	()	()	()	()
ICP cmH20	()	()	()	()	()	()	()	()	()
UAP mmHg	()	()	()	()	()	()	()	()	()
BAP mmHg	()	()	()	()	()	()	()	()	()
FAP mmHg	()	()	()	()	()	()	()	()	()
UVP mmHg	()	()	()	()	()	()	()	()	()
LV mmHg	()	()	()	()	()	()	()	()	()
Ao mmHg	()	()	()	()	()	()	()	()	()
T1/T2 ℃	36.6/36.5	36.6/36.4	36.6/36.4	36.6/36.5	36.6/36.4	36.6/36.4	36.6/36.5	36.6/36.4	36.6/36.4
CO2(Et/Fi) mmHg	38/2	38/2	38/2	38/2	38/2	38/2	38/2	38/2	38/2
02(Et/Fi) %	18/20	18/20	18/20	18/20	18/20	18/20	18/20	18/20	18/20
N20(Et/Fi) %	45/50	45/50	45/50	45/50	45/50	45/50	45/50	45/50	45/50



