

E65 ECG Machine Quick Reference Sheet

| Quick | Exam: | 66 | Power > Connect leads |
|---------------|---------------|---------|---|
| Shortcut Keys | | ÷ † - | Age Range |
| 10mm/mV | Waveform Gain | (†† | Gender |
| | | () | Mode Select > Auto |
| 25mm/s | Paper Speed | 0 | Paper Feed if needed |
| 100Hz | EMG Filter | | Record |
| | 4 14- | | Deserves |
| Patien | t Input: | | Power > Connect leads |
| Shortc | ut Keys | Patient | Patient Button LCD < Enter Data, press "OK" > |
| 10mm/mV | Waveform Gain | | Mode Select > Auto |
| 25mm/s | Paper Speed | 0 | Paper Feed if needed |
| 100Hz | EMG Filter | (2/0) | Record |

*Check paper and electrode supply daily