

12 lead ECG placement



- V₁** 4th intercostal space to the right of the sternum
- V₂** 4th intercostal space to the left of the sternum
- V₃** Directly between the leads V₂ and V₄
- V₄** 5th intercostal space at midclavicular line
- V₅** Level with V₄ at left anterior axillary line
- V₆** Level with V₅ at midaxillary line
(directly under the midpoint of the armpit)



Place limb leads on soft tissue surfaces rather than the bone according to this diagram:

- RA** Upper Right Arm
- LA** Upper Left Arm
- RL** Lower Right Leg
- LL** Lower Left Leg

Troubleshooting:

No Signal on Lead:

Check Electrode adhesion
Switch with next lead to eliminate electrode
If no signal evident after this, replace cableset.