





MAYO TROLLEY

Registration code: SKM-B





SAIKANG/Medical

Healthier your life

Product data sheet

SKH038-2

FEATURES AND HIGHLIGHTS

Technical parameters

 Size (LxWxH) 	500×350×(870-1300) mm
Tray size (LxWxH)	350×500×48 mm

Technical configuration

✓ Castor	4 pcs
✓ Tray	1 pc



Material

 $\Phi25mm,~~\Phi16mm$ stainless steel round pipe, the whole frame used high quality S.S. with the Mirror surface polishing.



Craftsmanship

Imported intelligent robot is used for welding forming, high precision, no blind spots, smooth surface, no burrs, beautiful and strong, high bearing capacity.



Tray

Size of stainless steel tray: L350*W500*H48mm.



Casto

a inch covered castors (diagonal brake), flexible and light to drive, smooth steering, wear-resistant, two with brake, can be easily braked in any state.



Product video (scan / press QR code)









Research



Manufacturing





Maintenance

Service



Saikang Medical Technology, Inc. 17800 Castleton St Ste 180, City of Industry, 91748, California. Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd. Dept. 111, 196 High Road Wood Green, N22 8HH, London. Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC 454080 Chelyabinsk. Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd. Viraj Complex Unit 3, Mombasa Road, P.O. Box 14262-00800, Nairobi. Kenya. Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC X22-S03 England Cluster, International City, Dubai. United Arab Emirates. Tel.: +20 128 8962497



Jiangsu Saikang Medical Equipment Co., Ltd. STOCK CODE: 870098

- No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
- □ +86 18021231901 · 24/7
- (the www.saikangmedical.com)





