



Designed with the latest in technical innovation, our feature-rich products deliver healthcare professionals a complete range of alternating pressure mattress solutions to suit long term and acute care environments.



- Advanced featured surface for high risk patients
- Optimal patient comfort and safety
- Easy set-up and automated therapy
- Comprehensive error and service indicators
- Cost effective solution



### **AutoCair**

Automatic weight sensing technology that delivers fail safe patient comfort levels. No more concerns over correct pressure level settings. Manual adjustments can still be made.



Cyclic



## Quicker Release CPR access

Easier quick twist mechanism located at the head of the mattress. Makes it easier to rapidly deflate the mattress for emergency CPR requirements.



**Quick Twist** 



## **Fowler Boost**

Adjusts mattress pressure to provide increased support under the sacral area when backrest is raised.



Manual



### SafetyCair

The ultimate protection from 'bottoming out', alternating air cells contain a lower 'safety cell' which remains fully inflated to provide greater surface stability and patient safety.




Air Cell



# King Single

Wider surface for therapy and positioning.







**Orange Medical Limited** PO Box 2070

Raumati Beach, Paraparaumu 5032

Phone: 0800 990024

Email: sales@orangenz.co.nz

www.orangenz.co.nz