

BVT01 | Bladder Volume Tester



Features

- 5.7"Color LED touch screen
- More accurate measurement, scan time<4s
- Easily review 12 frames of scan images
- Built-in thermal printer for easy printing
- USB port support massive data storage and software update
- Working time with Li-battery about 7 hours
- Multi-language customized(English,French,Russian,Spanish,Portuguese,Norway)

KAIXIN

XUZHOU KAIXIN ELECTRONIC INSTRUMENT CO.,LTD.

BVT01 | Bladder Volume Tester


Clinical Applications & Effects

Medical Services

- Urology
- Obstetrics-Gynecology
- Geriatrics
- Pediatrics
- Rehabilitation
- Surgical
- Operation Room
- Recovery Room
- Emergency
- Intensive Care
- Extended Care

Effects

- Diagnose urinary retention and evaluate many common urological conditions
- Prevent unnecessary catheterisation
- Reduce rates of urinary tract infection
- Monitor post-operative recovery
- Help caregivers manage and treat incontinence
- Reduce costs and save staff time



Technical Specifications

Description	Performance index	
Probe frequency, MHz	2.5	
Volume calculation reference plane	12plane , interval 15°	
Volume rendering resolution, mL	1mL	
Measuring range	20mL ~ 999mL	
Measuring accuracy, mL	20mL ~ 99mL	deviation±15mL
	100mL ~ 999mL	deviation±15%

Standard configuration

- | | |
|---------------|------|
| Main unit | 1pc |
| 2.5MHz probe | 1pc |
| Power adapter | 1pc |
| Battery | 1pc |
| Touch pen | 2pcs |

Optional configuration

- Battery charger
- Foot switch
- Trolley
- Verification cup



Tel: +86-516-87732932 / 87733758

Fax: +86-516-87732932 / 87792848

Http://www.kxele.com

E-mail: info@kxele.com / service@kxele.com

Zip Code: 221004

ADD: Kaixin Mansion, C-01, Economic Development Zone, Xuzhou, Jiangsu, China.



KAIXIN
XUZHOU KAIXIN ELECTRONIC INSTRUMENT CO.,LTD.